

**PP/SMS Isolation Gown (Classic Style)**

**Test Reports**



**AAMI Level 1**



**AAMI Level 2**



**AAMI Level 3**



**Meets AAMI Level 1, 2 & 3**

The gowns meet the Association for the Advancement of Medical Instrumentation (AAMI) PB70 standard, "Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities.

\*Light weight Gown---Level 1 : Protection against minimal levels of fluids

\*Medium weight Gown---Level 2 : Protection against low levels of fluids

\*High weight Gown---Level 3 : Protection against high levels of fluids



**10 pieces/bag**



**100 pieces/carton**  
Carton Size:52\*30\*34cm



**Warehouse**

Loading Quantity (25gsm SMS)		
1*20GP	1*40GP	1*40HQ
528 Cartons (52,800pcs)	1,056 Cartons(105,600pcs)	1,225Cartons(122,500pcs)

