



Fluid Repellent

Abrasion-Resistant

Low-Linting

Aseptic Folding

Ultrasonic Welding or Three thread over-lock Serged seam

Knitted Cuff

Latex free

Code	Product Name	Specification	Packaging
SGS10-35	Standard Gown	35gsm, Non sterile	5pcs/bag, 40pcs/case
SGS10-40	Standard Gown	40gsm, Non sterile	5pcs/bag, 40pcs/case
SGS10-45	Standard Gown	45gsm, Non sterile	5pcs/bag, 40pcs/case
SGS10-50	Standard Gown	50gsm, Non sterile	5pcs/bag, 40pcs/case
SGS11-35	Standard Gown	35gsm, sterile	1pc/sterile pouch, 25pcs/case
SGS11-40	Standard Gown	40gsm, sterile	1pc/sterile pouch, 25pcs/case
SGS11-45	Standard Gown	45gsm, sterile	1pc/sterile pouch, 25pcs/case
SGS11-50	Standard Gown	50gsm, sterile	1pc/sterile pouch, 25pcs/case

Note: Sterile gowns could include towel and wrapper in the package.